



MEDICATECH USA®

MasteRad® **MX Series** *Medical*



High Quality X-Ray Solutions
With a very Versatile Design

MasteRad™ MX30

Features

Fully integrated with all the features necessary for a all-in-one radiography system. Since your budget is always on our mind, MediatechUSA brings together a fully integrated system with smooth functions, user-friendly design and durable quality.



CONVENTIONAL
RADIOGRAPHY



FLOATING
TABLE



FLOOR WALL
TUBE RAIL
STAND



Ideal choice for a complete range of radiographic procedures and eliminates the need for high-cost room modifications. Enhance your workflow with a compact and effective design.

System Configuration & Specification

| | | |
|---------------------|--------------------------|-----------------------------|
| Generator | Input Power | 32 kW Single Phase |
| | kV Range | 40 -125 kV |
| | Output Frequency | 10 - 400mA |
| | mAs Range | 0.1 - 100 mA |
| | Tube Current Range | 10mA - 630mA 10mA - 800mA |
| | Exposure Time Range | 1 ms -6.4 s |
| | Console | Membrane Console |
| | A P R | Anatomical Programming |
| Technique Selection | kV/AEC, kV/mAs, kV/ma/ms | |

| | | |
|------|---------------------|----------------------|
| Tube | Type | Rotating Anode Type |
| | Focal Spot | Large 2mm, Small 1mm |
| | Operating Voltage | 40 - 150kV |
| | Anode Heat Capacity | 140kHu |
| | Target Angle | 12° |

| | | |
|------------|---------------------------|----------------------|
| Collimator | Power Supply | 1Phase 24VAC 50/60Hz |
| | Related Tube Voltage | 125 kVp |
| | Filters | 1.2mm Al |
| | Intensity of Illumination | LED over 160lux |


| | |
|------|---------------------|
| Grid | 103L 10:1 100/180cm |
| HVC | 10m |

| | | |
|------------|---------------------|------------------------------------|
| Tube Stand | Heights | 2250mm / 88.5 Inches |
| | Floor Rail Width | 3090mm / 121.5 Inches |
| | Vertical Movement | 400 - 1800mm / 15.75 - 70.5 Inches |
| | Horizontal Movement | 2433mm / 95.5 Inches |
| | Tube Stand Rotation | +/- 90° |
| | Tube Rotation | +/- 135° |
| | Telescopic Movement | 230mm / 9 Inches |

| | | |
|------------|--------------------|-------------------------------|
| Wall Stand | Wall Stand Heights | 2060mm / 81.25 Inches |
| | Vertical Movement | 355 - 1700mm / 14 - 70 Inches |

| | | |
|-------|-----------------|--|
| Table | Weight Capacity | 300kg |
| | Table Size | 2000x800x670mm / 78x31.5x26 Inches |
| | Movement | Transverse : 250mm / 9.75 Inches Longitudinal : 720mm / 28 Inches |
| | Top Filtration | 1.2mm / .04 Inches |

* All specifications can change without prior notice.



VOYANCE®

By Mediatech USA®

An industry leading complete software suite that delivers a secure and cost-effective way to acquire, view and share DICOM Images. Voyance® is one of the easiest, most user friendly interfaces. Its modern design & attention to the user experience along with its features & functions make it a great tool for every user.

Instant installation and auto-configuration

Easy, user-friendly Modern user interface

Interface consistency across acquisition and viewing platforms

Optimum Workflow designed for efficiency

ClairVoyance® Image Processing Algorithm

Supports a wide range of detectors and generators

X-RAY Generator integration (CPI & DRGEM)

Auto Stitching and Auto Cropping

Side by side comparisons

Fully DICOM 3.0 Compatible: MWL, MPPS, and all DICOM ops.

Reject Analysis





MEDICATECH USA®



T. +1 800-817-5030

T. +1 949-679-2881

F. +1 949-679-2882



www.medicattechusa.com



info@medicattechusa.com